



Head & Neck Cancer Rehab

Trismus, Radiation Fibrosis, Dysphagia, Xerostomia

Cleft * Craniofacial Disorders * Burns




Presurgical Positioning, Microstomia

Plastics * TMJD * Functional Rehab

Tissue Support, Post-Surgical, Jaw Ankylosis

Contents:

OraStretch® Press Jaw Rehab	2
OraStretch® Jaw Motion Scales	3
OraMPA™ Microstomia Appliance	3
DynaClose® Tissue System	3
OraPlace™ Facial Contracture Tools	4
MBI™ TMJ Imaging Prop	4
DynaCleft® PSIO System	5
nFlare™ Nasal Conformers	6
nPLANE™ Nasal Molding	6
Tongueometer™ Biofeedback Device	7
NutriSqueeze™ Feeding Bottles	7



OraStretch® Press Jaw Motion Rehab Device

Trismus * RT Fibrosis * TMJD * Stroke * Ankylosis * Burns * Surgical Rehab



Increase Oral Opening!

The **OraStretch®** Press jaw rehab system is a handheld passive motion device for the prevention and treatment of trismus and joint dysfunction.

We Confirm Benefits & Bill Insurance!


Medicare * BCBS * Aetna * Humana * TriCare & More

Features & Benefits:

- Wider opening for improved trismus prevention & recovery to full ROM.
- More versions plus free size exchanges*, for improved fit.
- Smaller size and weight. Easier to hold and use.
- Limited Lifetime Warranty*.



OraStretch User Kit



Dynamic Stretching with the Power Bands

OraStretch® Press Versions



Versions:	Population
Standard - CR1000	Adults; Children over 10 y.o.
Anterior - CR1003	For patients with access limitations.
Pediatric - CR1005	Children under 10 y.o.
Optional Additions:	Details:
Extended Biteplate	Clip-on biteplate extends 1/4". Severe over/under bites. Repeated pad pulling.
Edentulous Pads	3/8" Soft, thick pads for missing teeth, dentures. Cut to fit.
Power Bands	Provide Dynamic Stretching.

OraStretch® Scales - MIO, ROM, Lateral, Deviation



Unique Scales for your Patients

CR1085 - OraStretch MIO Scale - Paper Rotational - Industry standard with added micro-ROM and protrusion scales.

CR1087 - OraStretch Plastic MIO Scale - Single-patient reusable, telemedicine

CR1083 - OraStretch Deflection Scale - NEW! Clear middle to measure jaw deflection and deviation.

CR1080 - OraStretch ROM Scale - Paper Prong - Easily measure limited ROM

OraMPA™ Microstomia Appliance

LLPS Static, Dynamic & Manual Stretching

The OraMPA™ microstomia device offers static, dynamic and manual stretching of the lips and mouth for the prevention and treatment of microstomia. The device can also be used for orofacial muscle training. Standard, infant and XL prevention sizing.




Benefits:

- Easy-to-Use, Accessible Appliance
- Increased Lip and Mouth Opening ROM
- Multiple Therapy Options
- Low Load Prolonged Stretching (LLPS)


Microstomia Treatment for:

Orofacial Burns & Contractures, Scleroderma, Radiation Fibrosis, Post-Surgical Scarring, Lip Hypomobility, Craniofacial Disorders



DynaClose® Wound Closure and Tissue Expansion

Non-Invasive for Grafting or Wound/Excision Closure



Simple, non-invasive skin expansion for grafting, skin defect excision and grafting. DynaClose Strips work dynamically, moving as the body moves to provide a consistent appositional force that grows more skin quickly.

- Pre-op Stretch of the skin lesion area
- Expand Graft Sites
- Dynamic Wound Closure

Indications: Skin Defects, Retracted/Dehiscenced Wounds, Burns, Graft Expansion

DynaClose® Strips:

DWC24 - Small - 0.5" x 4.5" 15 pks of 4

DWC38 - Large - 0.75" x 8.5" 15 pks of 4



OraPlace™ Contracture Treatment and Prevention


Maintain Graft Size * Acute Microstomia Prevention * Increased ROM

The OraPlace™ Tools prevent and treat contractures with low load prolonged stretching (LLPS) of the face and mouth for improved range-of-motion (ROM), lip closure, and facial animation. Use for immediate intervention after burns or grafting, including with sedated and ventilated patients. Treat contractures to lengthen facial tissues for increased pliability and ROM.

**Facial Burns * Scarring * Radiation Fibrosis
Scleroderma * Cleft Lip * Upper Lip Tie Release**



Targeted Shapes in Graded Sizes



Bar

Target: Under philtrum and the mentolabial junction

Ball

Target: Deep pocket of nasal labial folds & oral commissures


Boomerang

Target: Cheeks and contours along the jaw line and upper molars

Benefits:

- Maintain and Support Full Graft Size
- Acute Intervention including on Sedated and Ventilated Patients
- Increase End-Range of Movement
- Improve ROM & Lip Closure
- Targeted Care for Facial Reanimation

Low Load Prolonged Stretching (LLPS)



MBI™ McCain-Burris Imaging Prop



Clear, easy, consistent open-bite MRI images

The MBI™ Imaging Prop simplifies & clarifies an open-mouth MRI of the TMJ. The MBI Prop offers three clinically relevant ROM options of 20mm, 30mm, and 40mm. The prop shows an open joint in rotation or translation, while reducing the opportunities to overstretch and injure the patient. The MBI Prop makes it easier for imaging technicians to produce valid usable images with simple and more accurate notation.

DynaCleft® Pre-Surgical System for Cleft Lip & Palate



Evidence-based treatment for babies with clefts

The DynaCleft® System gently guides the facial tissues to improve nasal symmetry and lip closure for simplified surgeries with improved cosmetic & functional outcomes.

Available in 7-packs or the New Baby Kits, which include supplies for up to the initial lip closure surgery. Use either option independently, together, or with an intra-oral plate. The DynaCleft system is an excellent alternative solution for PSIO when the patient does not tolerate NAM.

Learn more at: CranioRehab.com/DynaCleft



DynaCleft Positioning Strips

The strips gently apply to the cheeks to stretch & pull the cleft lip closer together for easier surgery. Unilateral or bilateral.

- Controlled force re-positions prolabium, premaxilla and alveolar segments.
- Add to NAM or intra-oral plate.



DynaCleft Nasal Elevators

The nasal clip supports the nose from a strap gently adhered to the forehead. Lifts and holds the nose up for balanced and even growth.

- Improves nasal symmetry and breathing.
- Lengthens and straightens the columella.
- Reduces need & extent of primary nasal surgery.

DynaCleft 7-Packs

DCX10 - Unilateral Strips **DCX20** - Bilateral Strips
DCNE-01 - Nasal Elevators

DynaCleft New Baby Kits

8 Strip + 6 Nasal Elevator packs for up to 12 weeks

DC1010 - Unilateral + Nasal Elevators
DC1015 - Unilateral Strips Only
DC1020 - Bilateral + Nasal Elevators
DC1025 - Bilateral Strips Only

Baby and Parent Benefits

- Pre-surgical positioning optimizes cosmetic & functional outcomes.
- Easy at-home treatment minimizes clinical visits - No doctor fittings necessary.
- Does not interfere with feeding.
- Non-invasive treatment maximizes comfort by stretching with the baby's mouth.
- Significant cost savings compared to NAM.



nFlare™ Nasal Stent & Conformer



Shape the Nose and Maintain the Airway

The nFlare™ Nasal Conformers help position and shape the nostrils and septum while maintaining an open airway. The nFlare conformers can be placed pre-surgically to shape the nostril & septum or post-surgically to stent and shape the tissues, support healing, and minimize scarring.

Post-Surgical Support & Shaping

- Nasal Reconstruction
- Rhinoplasties
- Cleft Surgery
- Septoplasties

Benefits:

- 15 Sizes for all ages and nasal sizes
- Symmetrical, Asymmetrical & Custom Shapes
- Shapes the Tissues and Maintains Opening
- Supports Healing and Minimizes Scarring
- Free Sizer for accurate selection



Asymmetric & Custom Shapes



nPLANE™ Rhinoplasty Appliance System

Pre- & Post-Surgical Lip and Nose Approximation

The nPlane device was designed to correct the patient's caudal nasal septal deviation, improve the nostril circumference and columellar length, elevate the nasal tip, and approximate the soft and hard tissue borders of the cleft defect.



Shape the Nose & Maintain the Airway

- 4 Sizes for Sequential Activation
- Improve & Maintain Nostril Shape
- Approximate Alveolar Segments
- Elongation of Columella
- Prevent Tissue Collapse
- Controlled Nasal Tip Projection

Tongueometer® Biofeedback Device

Tongue Strength and Endurance




Biofeedback at the tip of your tongue!

The **Tongueometer** device can be used to measure and increase tongue strength and endurance for patients with dysphagia.

The **Tongueometer** system works with a free downloadable Apple or Android app. The app features four modules that assess and treat tongue strength and endurance.

Features & Benefits:

- AFFORDABLE for use at home, on the go, or in the clinic
- Settings can be easily modified by the healthcare professional or user
- Stores exercise data within the app
- Easy sharing with clinicians
- Displays visual biofeedback



Assessment:

- Measure maximum tongue strength
- Measure patient's tongue endurance

Treatment:


- Track repetitions of tongue exercises
- Record isometric endurance exercises

Evaluation:

- Free downloadable excel template

NutriSqueeze™ Special Diet Bottles

Liquids, Thins, and Thicks! Simple with More Volume.



NutriSqueeze™ bottles allow patients to comfortably maintain adequate nutrition and hydration with improved access to specialty diets, easier feeding and medication delivery. We offer 4 straw types for your patients' needs. BPA, Latex, & DEHP Free

NutriSqueeze Bottles

NutriSqueeze™ Cleft Bottle - 4oz: CR6014 - Silicone Straw with volume info

NutriSqueeze™ Liquid Bottle - 16oz: CR6016, 8oz: CR6018 - Hard curved straw

NutriSqueeze™ Thins Bottle - 8oz - CR6023: 3/16" Silicone Straw CR6024: 1/4" PVC Straw

NutriSqueeze™ Thicks Bottle - 16oz: CR6026, 8oz: CR6028 - 3/8" Soft thick straw





See product IFU or call us for warnings, cautions, indications, contraindications, and instructions of use for products shown.

Any claims are general benefits and results may vary.

OraStretch®, TheraPacer™, OraMPA™, OraPlace™, MBI™, nFlare™, NPlane™, Salvate™, Hiberna™ & NutriSqueeze™ are trademarks of CranioMandibular Rehab, Inc. All other trademarks are property of their respective owners. ©CranioMandibular Rehab, Inc.



ProCMR2411



CranioRehab International Products:

Trismus and Orofacial Fibrosis

OraStretch® Press Jaw Motion Rehab System & Scales

Jaw exercise and stretching to treat and prevent trismus, joint dysfunction and improved surgical rehab.

OraPlace™ Contracture Tools - Localized orofacial scarring and fibrosis treatment.

OraMPA™ Microstomia Appliance - Prevent and treat microstomia and lip contractures.

MBI™ McCain-Burris Imaging Prop - Improve TMJ open-bite MRI images.

Nasal & Cleft Management

DynaCleft® Positioning Strips & Nasal Elevators

Gently guide & support facial tissues for improved cleft lip and palate surgery outcomes.

nFlare™ Nasal Stent Conformer - Stent, shape, maintain and expand the nostrils.

nPlane™ RAS Nasal Molding - Expand and shape the nostrils while maintaining the airway.

Dysphagia & Xerostomia

Tongueometer™ Biofeedback Device - Measure and increase tongue strength & endurance.

Pre- and Post-Surgical Care

NutriSqueeze™ Bottles - Specialty feeding bottles for liquids, thins, and thicks.

DynaClose® Tissue Expansion & Wound Closure

Non-invasive skin expansion, wound closure, dehiscence prevention, and suture support.

